

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

TGM4 Anti-Human Antibody (Clone AT1H1)

单克隆抗体

产品货号: MA01332

第三版

免疫原:Recombinant human TGM4 (421-565aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:AT1H1

描述:TGM4 (Transglutaminase 4; also TGP and Prostate-specific transglutaminase) is a 77kDa member of the papain-like family of transglutaminases. TGM4 is a unique member of the family due largely to its distribution pattern in the body in that it is primarily found in the prostate gland and existed only at low levels in other tissues types.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF