

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Tight Junction-Associated Protein (7H6) Mouse IgM MoAb

单克隆抗体

产品货号: MA01336

第三版

免疫原:A bile canaliculus-rich fraction prepared from rat livers

亚型:IgM

克隆:7H6

描述:Tight junctions function to maintain cellular polarity and permeability barriers in epithelial and endothelial cells. Recently, several component molecules of tight junction have been identified.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF