

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

TLR7 Anti-Human Antibody (Clone 4F4)

单克隆抗体

产品货号: MA01337

第三版

免疫原: Recombinant human TLR7 (451-500aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: 4F4

描述: TLR7 is one of toll-like receptors which are type I transmembrane proteins involved in innate immunity by recognizing microbial conserved structures. These proteins are characterized by multiple copies of leucine rich repeats in the extracellular domain and IL-1 receptor motif in the cytoplasmic domain, and act through adaptor molecules such as MyD88 and TIRAP to activate various kinases and transcription factors.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF