
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

TMEFF2 Anti-Human Antibody (Clone J4B6)

单克隆抗体

产品货号: MA01338

第三版

免疫原:Recombinant human TMEFF2 (41-150aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:J4B6

描述:TMEFF2(transmembrane protein with EGF-like and two follistatin-like domains 2), a gene encoding a plasma membrane protein with two follistatin-like domains and one epidermal growth factor – like domain, had limited normal tissue distribution and was highly overexpressed in prostate cancer.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF