

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Tob (4B1) Mouse IgG MoAb

单克隆抗体

产品货号: MA01341

第三版

免疫原:Recombinant Tob (E. coli)

亚型:IgG2a

克隆:4B1

描述:A cDNA for a novel protein termed Tob (Transducer of ErbB-2) that interacts with the c-erbB-2 gene product p185erbB2 was molecularly cloned. Nucleotide sequencing reveals that the Tob protein is a 45 kDa protein that does not contain either SH2 or SH3 domain but is homologous to the previously characterized anti-proliferative gene product BTG-1 at its amino-terminal half.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF