



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Tollip Anti-Human Antibody (Clone K1B17)

单克隆抗体

产品货号: MA01343

第三版

**免疫原:** Recombinant human Tollip (61-274aa) purified from E. coli

**亚型:** Mouse IgG2a heavy chain and kappa light chain

**克隆:** k1B17

**描述:** Toll interacting protein, also known as Tollip, is an inhibitory adaptor protein within the signaling pathway of IL-1 and Toll-like receptors (TLR). The Toll-like receptors are essential transmembrane signaling receptors of the innate immune system that alert the host to the presence of a microbial invader.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF