
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

TPD52L1 Anti-Human Antibody (Clone d1C5)

单克隆抗体

产品货号: MA01347

第三版

免疫原:Recombinant human TPD52L1 (1-131aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:d1C5

描述:TPD52L1(hD53,hD53L1) is a member of D52 gene family and overexpressed in human breast carcinoma. This protein contains coiled-coil motif, which is responsible for homo- and heteromeric interaction, and may be a regulator of cell proliferation and calcium-mediated signal transduction.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF