



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

TRA1 Anti-Human Antibody (Clone 2H3)

单克隆抗体

产品货号: MA01349

第三版

免疫原: Recombinant human TRA1 (676-803aa) purified from E. coli

亚型: Mouse IgG2a heavy chain and kappa light chain

克隆: 2H3

描述: Tumor rejection antigen 1 (TRA1) is expressed on the surface of human teratocarcinoma stem cells (EC), human embryonic germ cells (EG) and human embryonic stem cells (ES). TRA1 is a marker protein of these stem cells.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF