



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Trem2 Anti-Human Antibody (Clone 2B5)

单克隆抗体

产品货号: MA01350

第三版

免疫原: Recombinant human TREM2 (19-161aa) purified from E. coli

亚型: Mouse IgG3 heavy chain and kappa light chain

克隆: 2B5

描述: The triggering receptor expressed on myeloid cells-2 (TREM2) is a membrane protein that interacts with the TYROBP protein. The TREM2 and TYROBP proteins form a complex that transmits chemical signals to activate the cell.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF