

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

TTR Anti-Human Antibody (AT7B2)

单克隆抗体

产品货号: MA01353

第三版

免疫原: Recombinant human TTR (21-147aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: AT7B2

描述: TTR encodes transthyretin, a circulating and abundant tetramer that functions as a carrier protein of thyroid hormone and retinol binding protein. Transthyretin (TTR) is mainly produced in the liver, but also in the brain and eye.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF