

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Tubb2B Anti-Human Antibody (At5B3)

单克隆抗体

产品货号: MA01354

第三版

**免疫原:** Recombinant human TUBB2b (1-445aa) purified from E. coli

**亚型:** Mouse IgG2a heavy chain and kappa light chain

**克隆:** AT5B3

**描述:** The protein encoded by this gene is a beta isoform of tubulin, which binds GTP and is a major component of microtubules. Multiple beta Tubulin isoforms have been characterized and are expressed in mammalian tissues.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF