

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

TUBG1 Antibody (Clone AT1C12)

单克隆抗体

产品货号: MA01355

第三版

免疫原:Recombinant human TUBG1 (1-451aa) purified from E.coli.

亚型:Mouse IgG2b, k

克隆:AT1C12

描述:Tubulin is the major constituent of microtubules. Tubulin, gamma 1 is a protein in humans that is encoded by the TUBG1 gene. This gene encodes a member of the tubulin superfamily. The encoded protein localizes to the centrosome where it binds to microtubules as part of a complex referred to as the gamma-tubulin ring complex. The protein mediates microtubule nucleation and is required for microtubule formation and progression of the cell cycle.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF