



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

-Tubulin Anti-Human Antibody (Clone AT5B2)

单克隆抗体

产品货号: MA01356

第三版

免疫原:Recombinant human beta-Tubulin (1-445aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:AT5B2.

描述:Tubulin is a major cytoskeleton component that has five distinct forms. Multiple beta Tubulin isoforms have been characterized and are expressed in mammalian tissues. Tubulin is the major constituent of microtubules.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF