



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Tubulin Beta 3 Antibody (Clone AT2F5)

单克隆抗体

产品货号: MA01357

第三版

**免疫原:**Recombinant human TUBB3 (1-450aa) purified from E. coli

**亚型:**Mouse IgG1, k

**克隆:**AT2F5

**描述:**Tubulin Beta 3 is 450 amino acid protein. Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha-chain. TUBB3 plays a critical role in proper axon guidance and maintenance. Its expression is primarily restricted to central and peripheral nervous system.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在-20 °C或者-70 °C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF