



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

TYMS Anti-Human Antibody (Clone AT21S5)

单克隆抗体

产品货号: MA01359

第三版

**免疫原:**Recombinant TYMS (1-313aa) purified from E. coli

**亚型:**Mouse IgG2b heavy chain and kappa light chain

**克隆:**AT1S5

**描述:**The methylation of deoxyuridine monophosphate (dUMP) to deoxythymidine monophosphate (dMTP) is an essential step for the formation of thymine nucleotides. This process is catalyzed by thymidylate synthase (TS or TYMS).

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF