

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

UBE2L6 Anti-Human Antibody (Clone k1H3)

单克隆抗体

产品货号: MA01360

第三版

免疫原:Recombinant human UBE2L6 (1-152aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:k1H3

描述:Ubiquitin-conjugating enzyme E2L6 (UBE2L6), also known as UbcH8, is a member of the E2 ubiquitin-conjugating enzyme family. The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在-20 °C或者-70 °C. 避免反复冻融.

应用:ELISA, WB, ICC/IF