
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

UBE2S Anti-Human Antibody (Clone AT2C12)

单克隆抗体

产品货号: MA01361

第三版

免疫原:Recombinant human UBE2S (1-222aa) purified from E. coli

亚型:Mouse IgG2a heavy chain and kappa light chain

克隆:AT2C12

描述:UBE2S, also known as ubiquitin-conjugating enzyme E2S, is a member of the ubiquitinconjugating enzyme family. This protein is able to form a thiol ester linkage with ubiquitin in an ubiquitin activating enzyme-dependent manner, a characteristic property of ubiquitin carrier proteins.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF