

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

UNG Anti-Human Antibody (Clone k1C12)

单克隆抗体

产品货号: MA01363

第三版

免疫原:Recombinant human UNG (1-313aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:k1C12

描述:The human UNG gene encodes both mitochondrial (UNG1) and nuclear (UNG2) forms of uracil- DNA glycosylase (UNG). These forms are generated from transcription from alternative promoters, promoter A and promoter B respectively, and the subsequent use of alternative splicing.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF