
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Vaspin Anti-Human Antibody (Clone 1C4)

单克隆抗体

产品货号: MA01364

第三版

免疫原:Recombinant human Vaspin (21-414aa) purified from E. coli

亚型:Mouse IgG1 heavy chain and kappa light chain

克隆:1C4

描述:Vaspin (visceral adipose-specific SERPIN), a newly identified adipokine, which is a member of serine protease inhibitor family. Vaspin is also a unique insulin sensitizing adipocytokine in obesity.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF