



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

VEGF (16F1) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01365

第三版

免疫原:Recombinant Human VEGF (E. coli)

亚型:IgG2a

克隆:16F1

描述:Vascular Endothelial Cell Growth Factor (VEGF) is a homodimeric protein initially purified from media conditioned by normal bovine pituitary folliculo-stellate cells and secreted by a variety of vascularized tissues.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF