

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-VIP Receptor Rat IgG MoAb

单克隆抗体

产品货号: MA01373

第三版

免疫原: Human adenocarcinoma cells (HT29 cell line) originating from human colon.

亚型: IgG2c

克隆: 109-10

描述: This antibody can be stained in frozen section. The optimal dilution is X100, however, the dilution rate should be optimized by each laboratories.

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度: Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF