



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Visfatin Anti-Human Antibody (Clone 4D5)

单克隆抗体

产品货号: MA01374

第三版

**免疫原:** Recombinant human Visfatin (1-491aa) purified from E. coli

**亚型:** Mouse IgG1 heavy chain and kappa light chain

**克隆:** 4D5

**描述:** Visfatin is pre-B cell colony-enhancing factor 1 (PBEF), which is predominantly secreted from visceral adipose tissue both in humans and mice. This protein has also been reported to be a cytokine (PBEF) that promotes B cell maturation and inhibits neutrophil apoptosis, or a visceral-fat derived hormone that acts by binding and activating the insulin receptor.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF