



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Wnt3a Anti-Human Antibody (Clone 3A6)

单克隆抗体

产品货号: MA01375

第三版

**免疫原:** Recombinant human Wnt3a (19-352aa) purified from E. coli

**亚型:** Mouse IgG2a heavy chain and kappa light chain

**克隆:** 3A6

**描述:** Wingless-type MMTV integration site family, member 3A (Wnt3a) is a member of the Wnt family. The Wnt gene family is representative of a growing class of multifunctional factors that play a role in both tumorigenesis and patterning events required for normal differentiation and development.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF