

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

XPNPEP1 antibody (clone AT9C7)

单克隆抗体

产品货号: MA01377

第三版

**免疫原:**Recombinant human XPNPEP1 (1-623aa) purified from E. coli

**亚型:**IgG2b, k

**克隆:**Anti-human XPNPEP1 mAb, clone AT9C7, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human XPNPEP1 protein.

**描述:**XPNPEP1 also known as Xaa-Pro aminopeptidase 1 isoform 1, is encodes the cytosolic form of a metalloaminopeptidase that catalyzes the cleavage of the N-terminal amino acid adjacent to a proline residue. This product may play a role in degradation and maturation of tachykinins, neuropeptides, and peptide hormones.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF