
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

YB-1 (21A3) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01378

第三版

免疫原:Recombinant human YB-1

亚型:IgG1

克隆:21A3

描述:YB-1 (Y-box binding protein-1) contains a highly conserved cold shock domain and is a transcription factor that binds to a CCAAT box. It has been reported that P-glycoprotein is highly expressed in cisplatin-resistant cancer cell lines since the sequence of CCAAT box is located in the promoter region of the drug-resistance gene MDR1, and YB-1 is involved in transcription of P-glycoprotein, the product of the MDR1 gene.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

浓度:100 μg/mL

保存: +4 °C 保存 (1-2 周). 长期保存在-20 °C或者-70 °C. 避免反复冻融.

应用:ELISA, WB, ICC/IF