

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Hexokinase 1, 2 Anti-Human Antibody (Clone 1A7)

单克隆抗体

产品货号: MA01412

---

第三版

**免疫原:** Recombinant human Hexokinase 2 (1-917aa) purified from E. coli

**亚型:** Mouse IgG1 heavy chain and kappa light chain

**克隆:** Anti-human Hexokinase 1, 2 mAb, clone 1A7, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

**描述:** Hexokinase is the first enzyme in the glycolytic pathway, catalyzing the transfer of a phosphoryl group from ATP to glucose to form glucose-6-phosphate and ADP. In mammals, four distinct enzymes (types 1 to 4 hexokinases) have been identified.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF