

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-HIV-1 p24 (2C2) Mouse IgG MoAb

单克隆抗体

产品货号: MA01414

第三版

免疫原:Recombinant HIV-1 gag p24

亚型:IgG1

克隆:2C2

描述:HIV-1 is the retrovirus which is cause virus of the AIDS. HIV-1 p24 is a capsid protein constituting the core of HIV-1. HIV-1 p24 protein is detectable from the early days of HIV-1 infection. HIV-1 p24 has been measured as a marker of posttreatment pathologic condition to know quantity of virus in blood.

配方:Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN₃

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF