



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-HIV-1 p24 (NP-24) Mouse IgG MoAb

单克隆抗体

产品货号: MA01415

第三版

免疫原: Recombinant HIV-1 gag p24

亚型: IgG1

克隆: NP-24

描述: HIV-1 is the retrovirus which is cause virus of the AIDS. HIV-1 p24 is a capsid protein constituting the core of HIV-1. HIV-1 p24 protein is detectable from the early days of HIV-1 infection. HIV-1 p24 has been measured as a marker of posttreatment pathologic condition to know quantity of virus in blood.

配方: Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN<sub>3</sub>

浓度: Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF