



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Homer1 Anti-Human Antibody (Clone AT1F3)

单克隆抗体

产品货号: MA01418

第三版

免疫原: Recombinant human Homer1 (1-354aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: Anti-human Homer1 mAb, clone AT1F3, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述: Homer1 (Homer protein homolog 1) belongs to Homer family of proteins that regulate signal transduction through, and the trafficking of, glutamate receptors, as well as maintain and regulate extracellular glutamate levels in corticolimbic brain regions. The Homer1 is enriched at excitatory synapses and binds group 1 metabotropic glutamate receptors (mGluRs)

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF