

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Hsf1 Anti-Human Antibody (Clone J7F9)

单克隆抗体

产品货号: MA01419

第三版

免疫原: Recombinant human HSF1 (1-529aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: Anti-human HSF1 mAb, clone J7F9, is derived from hybridization of mouse FO myeloma cells with spleen cells from BALB/c mice

描述: Heat shock factor 1 (HSF1) is a major transactivator of heat shock proteins in response to environmental changes, and it is also involved in ogenesis, spermatogenesis, and placental development.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF