

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Hsp105 Antibody (Clone J1G12)

单克隆抗体

产品货号: MA01420

第三版

免疫原: Recombinant human Hsp105 (1-858aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: J1G12.

描述: Hsp105 is a mammalian member of the Hsp105/110 family, a diverged subgroup of the Hsp70 family. Hsp105 exists as two isoforms, alpha and beta. Hsp105 alpha associates with Hsp70/Hsc70 as complexes in vivo and regulates the chaperone activity of Hsp70/Hsc70 negatively in vitro and in vivo.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF