



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Hsp40 Anti-Human Antibody (Clone k1C7)

单克隆抗体

产品货号: MA01422

第三版

**免疫原:** Recombinant human Hsp40(1-340aa) purified from E. coli

**亚型:** Mouse IgG2a heavy chain and kappa light chain

**克隆:** k1C7

**描述:** Human Hsp40 (heat shock protein with molecular size of 40kDa) is one of the mammalian homologues of bacterial DnaJ heat shock protein and regulates cellular processes by aiding in the folding, transport and assembly. Hsp40 protein contains the J-domain which controls interaction with the ATPase domain of DnaK.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF