

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Hsp70 Anti-Human Antibody (Clone 4E7)

单克隆抗体

产品货号: MA01424

---

第三版

**免疫原:** Recombinant human Hsp70 (1-641aa) purified from E. coli

**亚型:** Mouse IgG1 heavy chain and kappa light chain

**克隆:** 4E7.

**描述:** Hsp70 is a human heat shock protein. Hsp70 is an important part of the cell's machinery for protein folding, and help to protect cells from stress. In most species, there are many proteins that belong to the hsp70 family.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF