



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Hsp90 Anti-Human Antibody (Clone 4F10)

单克隆抗体

产品货号: MA01425

第三版

免疫原: Recombinant human Hsp90 (1-732aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: 4F10.

描述: Hsp90 is a human heat shock protein. In response to adverse change in their environment, cell from all organisms increase the expression of a class of proteins referred to as heat shock or stress protein.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF