



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

HtrA2/Omi Anti-Human Antibody (Clone 1B3)

单克隆抗体

产品货号: MA01427

第三版

免疫原: Recombinant human HtrA2/Omi (134-458aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: 1B3.

描述: HtrA2/Omi is a mammalian serine protease at high temperatures and a chaperone at low temperature. The full-length HtrA2 is synthesized as a precursor protein and then targeted to the mitochondria where it is matured by the removal of N-terminal 133 residues. Recently, HtrA2 has known to contribute both to caspase-dependent and caspase-independent cell death

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF