

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

IDH1 R132H (HMab-1) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01431

---

第三版

**免疫原:**Synthetic peptide of a part of human IDH1 R132H

**亚型:**IgG1

**克隆:**HMab-1

**描述:**Specific mutations in the isocitrate dehydrogenase 1 gene IDH1 have been found in several brain tumors including astrocytoma, oligodendroglioma and glioblastoma multiforme, with mutations found in nearly all cases of secondary glioblastomas, but rarely in primary high-grade glioblastoma multiforme.

**配方:**Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN<sub>3</sub>

**浓度:**Please see datasheet.

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF