
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

IDH1 R132S (SMab-1) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01433

第三版

免疫原: Synthetic peptide of a part of human IDH1 R132S

亚型: IgG1

克隆: SMab-1

描述: Specific mutations in the isocitrate dehydrogenase 1 gene IDH1 have been found in several brain tumors including astrocytoma, oligodendroglioma and glioblastoma multiforme, with mutations found in nearly all cases of secondary glioblastomas, but rarely in primary high-grade glioblastoma multiforme.

配方: Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃

浓度: Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF