



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

IL-32 Anti-Human Antibody (Clone AT2F9)

单克隆抗体

产品货号: MA01437

第三版

免疫原: Recombinant human IL32 (1-131aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: AT2F9

描述: The gene encoding IL-32 (Interleukin-32) is located on human chromosome 16p13.3 and is organized into eight exons. IL32 is mainly expressed in natural killer cells, T cells, epithelial cells, and blood monocytes.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF