

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

IMP3 Antibody (Clone AT32E9)

单克隆抗体

产品货号: MA01440

第三版

**免疫原:** Recombinant human IMP3 (1-184aa) purified from E. coli

**亚型:** IgG1, k

**克隆:** Anti-human IMP3 mAb, clone AT32E9, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human IMP3 protein.

**描述:** IMP3 is the human homolog of the yeast Imp3 protein. The protein localizes to the nucleoli and interacts with the u3 snoRNP complex. The protein contains an S4 domain.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF