

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Intelectin-1 (3G9) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01442

第三版

免疫原:Recombinant human Intelectin-1

亚型:IgG2a

克隆:3G9

描述:Intelectin-1 is a secretory animal lectin that binds to galactofuranose. Human intelectin-1 is a glycoprotein mainly secreted from intestinal goblet cells.

配方:Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN<sub>3</sub>

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF