



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Interferon Regulatory Factor-3 Antibody (Clone 3F10)

单克隆抗体

产品货号: MA01444

第三版

免疫原: Recombinant human IRF-3(108-166aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: 3F10

描述: Members of the Interferon regulatory factor (IRF) family regulate gene expression critical to immune response, hemopoiesis, and proliferation. IRF-3 is a member of the IRF family, and is distinct from other family members. Its transcriptional activity is regulated solely by posttranslational modifications. It plays a crucial role in activation of innate immunity and inflammation in response to viral infection.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF