

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

IRF-5 Anti-Human Antibody (Clone 10T1)

单克隆抗体

产品货号: MA01445

第三版

免疫原: Recombinant human IRF-5 (25-235 aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: 10T1

描述: Interferon regulatory factor 5 (IRF-5) is a member of IRF family, which plays a critical role in the regulation of genes induced by viral infection, immunostimulation and cell growth regulation.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF