

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Irf7 Anti-Human Antibody (Clone 3D9)

单克隆抗体

产品货号: MA01446

第三版

免疫原:Recombinant human IRF-7 (1-150aa) purified from E. coli

亚型:Mouse IgG1 heavy chain and kappa light chain

克隆:3D9

描述:Interferon regulatory factor (IRF) 7 is a member of the interferon regulatory transcription factor family. IRF-7 plays an essential role in virus-activated transcription of type I interferon alpha genes. Inducible expression of IRF-7 is largely restricted to lymphoid tissue. Multiple IRF-7 transcript variants have been identified, although the functional consequences of these have not yet been established.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF