

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ISG15 Anti-Human Antibody (Clone 3E5)

单克隆抗体

产品货号: MA01447

第三版

免疫原: Recombinant human ISG15 (1-157aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: 3E5.

描述: ISG15, a 15 kDa protein encoded by an interferon stimulated gene, and its modification system play important roles in the innate immune response, regulation of interferon signaling, pregnancy, and several cancers. ISG15 is an ubiquitin-like protein containing two ubiquitin homology domains and becomes conjugated to a variety of proteins when cells are treated with type I interferon or lipopolysaccharide.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF