

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ITLN1 Antibody (Clone AT6C11)

单克隆抗体

产品货号: MA01448

第三版

免疫原: Recombinant human ITLN1 (17-313aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and light chain

克隆: Anti-human ITLN1 mAb, clone AT6C11, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human ITLN1 protein.

描述: ITLN1, also known as intelectin1 or omentin, has no effect on basal glucose uptake but enhances insulin-stimulated glucose uptake in adipocytes. ITLN1 may play a role in the defense system against microorganisms and specifically recognize carbohydrate chains of pathogens and bacterial components containing galactofuranosyl residues, in a calcium-dependent manner.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF