

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

JAK2 Antibody (Clone AT53B7)

单克隆抗体

产品货号: MA01449

第三版

**免疫原:** Recombinant human JAK2 (1014-1132aa) purified from E. coli

**亚型:** Mouse IgG2b, k

**克隆:** Anti-human JAK2 mAb, clone AT53B7, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human JAK2 protein.

**描述:** JAK2 is a protein tyrosine kinase involved in a specific subset of cytokine receptor signaling pathways. It has been found to be constitutively associated with the prolactin receptor and is required for responses to gamma interferon. Mice that do not express an active protein for this gene exhibit embryonic lethality associated with the absence of definitive erythropoiesis.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF