

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

KCTD15 Anti-Human Antibody (AT4C3)

单克隆抗体

产品货号: MA01450

第三版

免疫原:Recombinant human KCTD15 (1-234aa) purified from E. coli

亚型:Mouse IgG3 heavy chain and L light chain

克隆:AT4C3

描述:KCTD15 (Potassium channel tetramerisation domain containing 15, also known as BTB/POZ domain-containing protein KCTD15) is protein that in humans is encoded by the KCTD15 gene. KCTD15 is expressed at high level in brain and hypothalamus.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF