



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

KIR2DL1 Anti-Human Antibody (Clone 2F9)

单克隆抗体

产品货号: MA01460

第三版

**免疫原:** Recombinant human KIR2DL1 (23-223 aa) purified from E. coli

**亚型:** Mouse IgG2a heavy chain and kappa light chain

**克隆:** 2F9

**描述:** The KIR family consists of transmembrane glycoproteins of the Ig superfamily expressed on human NK cells and a subset of human T cells which they are involved in recognition of either MHC class I molecules or unknown ligand on target cells and inhibit cytotoxic activities.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF