



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

KIR2DL3 Anti-Human Antibody (Clone190IIC311)

单克隆抗体

产品货号: MA01461

第三版

**免疫原:**Recombinant human KIR2DL3(19-161aa) purified from E. coli

**亚型:**Mouse IgG2a heavy chains and light chain

**克隆:**190IIC311

**描述:**Reacts with extracellular domain of the KIR2DL3, a member of the killer cell immunoglobulin-like receptor(KIR) family. The KIR family consists of transmembrane glycoproteins of the Ig superfamily expressed on human NK cells and a subset of human T cells which they are involved in recognition of either MHC class I molecules or unknown ligand on target cells and inhibit cytotoxic activities.<sup>2</sup>

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF