

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

KIR2DS4 Anti-Human Antibody (Clone 5F2)

单克隆抗体

产品货号: MA01462

第三版

免疫原: Recombinant human KIR2DS4 purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: 5F2

描述: The KIR family consists of transmembrane glycoproteins of the Ig superfamily expressed on human NK cells and a subset of human T cells which they are involved in recognition of either MHC class I molecules or unknown ligand on target cells and inhibit cytotoxic activities. KIR2DS4 is an activating receptor of KIR family.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF